

Multi-animal Monitoring and Gating Systems for use with MR, PET, CT, SPECT and Optical

Configurations:

- MR
- Imaging systems other than MR
- 3 or 4 channels

Temperature control

Waveform & trend data acquisition



Monitoring

- ECG
- Respiration
- Temperature

Retrospective Gating

- ECG
- Respiration
- ECG & respiration

Multi-animal monitoring and gating systems have been designed to meet the physiological monitoring and gating needs for anesthetized mice, rats and larger animals in the MR, PET, CT, SPECT, Optical and laboratory environments. The Model 1200 is configured with a 3 or 4 channel ECG, respiration and temperature (ERT) data acquisition Module located near the imaging volume and near the animals and with a PC Interface (PCI) Module connected to a PC located near the operator console. The PC displays multiple waveforms, measured values, trends and gating pulses. Data acquisition is controlled by menu driven software from the PC.

ECG waveforms are measured for each animal using two or three leads with sub-dermal needle, gold disk surface or radio translucent surface electrodes. Each ECG waveform is processed to detect the R-wave, generate an ECG gate and to determine the heart rate.

Temperature is measured using either thermistor or fiber optic sensors. Fiber optic sensors offer several advantage but are more expensive than thermistor sensors. The sensor can be located in the rectum or on the animal's skin. Any animal temperature measurement can be used with the Heater System to control the temperature of the animals.

Respiration is obtained from a small pneumatic pillow sensor placed next to the abdomen or in MR from one of the ECG leads for each animal. Each respiration waveform is processed to detect inspiration, expiration, a respiration gate and respiration rate.

Auxiliary input channels allow the user to synchronize the animal's physiological measurements with data from the imaging system which allows retrospectively gated images.

Gating algorithm, which is user configured, allows gates to be generated for each animal from ECG, respiration or ECG and respiration. The user can control the start, stop and width of each animal's gate.


Power for the ERT Module is supplied from either a battery or an isolated power supply.

Specifications:

ECG	range: accuracy: input range: input Impedance: CMRR:	40 - 900 BPM $\pm 1\%$ -2.50 mV to 2.5mV >10 M Ω at 10 Hz 100 dB at 60 Hz
Temp	probe types tip diameter range accuracy	thermistor or fiber optic 1.0 or 3.0 mm 20 - 60 $^{\circ}\text{C}$ ± 0.26 $^{\circ}\text{C}$ 32 - 42 $^{\circ}\text{C}$
Resp	range accuracy sensor type	13 - 300 bpm 1 count pneumatic pillow or ECG lead
ERT Module	multiple channels power patient isolation	3 or 4 battery or power supply optical
PCI Module	aux gate inputs power	3 or 4 +12 VDC
Gating	R-wave to gate delay expiration gate width and delay	selectable - 0 ms to 600 ms selectable - 1 ms step size
Temp	heater control	fiber optic PWM

PC requirements:

Operating system:	Windows any version
Hardware:	>1 GHz processor Display $\geq 1024 \times 768$ USB port




SA Instruments, Inc.

65 Main Street, Stony Brook, NY 11790 www.i4sa.com

(631) 689-9408

FAX: (631) 689-9410

Model 1200 Multi-animal Monitor Display



cardiac gate:
red dot when
R-wave is
detected,
white dot when
gate is output

respiratory
gate: red bar
during
expiration,
white dot when
gate is output

ECG, respiration waveforms and temperature trend displayed for the selected channel. Heart, respiration rates and temperature displayed for all channels.



SA Instruments, Inc.
Stony Brook, NY